INSTRUCTIONS FOR USE
Gynecological Kit

【Product Name】 Gynecological Kit

【Product Model】

HH850401, HH850402, HH850403, HH850404, HH850405, HH850406, HH850407, HH850408, HH850409, HH850410, HH850411, HH850412, HH850413, HH850414, HH850415, HH850416, HH850417, HH850418, HH850419, HH850420, HH850421, HH850422, HH850423, HH850424, HH850425, HH850426, HH850427, HH850428, HH850429, HH850430, HH850431, HH850432, HH850433, HH850434, HH850435, HH850436, HH850437, HH850438, HH850439, HH850440, HH850441, HH850442, HH850443, HH850444, HH850445, HH850446, HH850447, HH850448, HH850449, HH850450, HH850451, HH850452, HH850453, HH850454, HH850455, HH850456, HH850457, HH850458, HH850459, HH850460, HH850461, HH850462, HH850463, HH850464, HH850465, HH850466, HH850467, HH850468, HH850469, HH850470, HH850471, HH850472

【Product structure】

The gynecological kit consists of vaginal dilator, sampling brush, sampling swab, cervical sampler, examination glove, medical drape and microscope slide (optional). The configurations of each model are different.

【Intended Use】

It is mainly used for gynecological examination.

【Contraindication】

1) Disabled within 2 days before and after menstruation.

2) pregnant women are banned.

3) Vaginal bleeding is prohibited.

【Steps to use】

1. Open the product pouch packaging.

2. Spread the medical drape on the examination bed.

3. The patient lies on the medical drape, takes off her pants and places her legs on the bracket of the examination bed.

4. Gently insert the vaginal dilator into the patient's vagina and then slowly open it.

5. Insert the head of the sampling swab into the vagina and rotate it to remove vaginal secretions.

6. Insert the head of the sampling brush or the head of the cervical sampler to the vaginal sampling area and rotate it 2 ～ 5 turns in one direction for sampling.

7. After sampling is completed, evenly spread the sample on the microscope slide or place the sampling head into the preservation solution for subsequent testing.

【Precaution】

1. Please check the pouch packaging and the product carefully before use. If either of them is damaged, it is strictly forbidden to use.

2. Please read the instruction manual carefully, pay attention to the expiry date of the product before use. Do not use it after the deadline.

3. This product is disposable. It should be strictly in accordance with the aseptic operation specifications. Please destroy it immediately after use. It is forbidden to reuse.

4. The disposal of wastes shall be carried out in accordance with national environmental protection laws and regulations.

【Definitions of Signs】

STERILE EO	Sterilized using ethylene oxide		Do not re-use
	Do not use if package is damaged		Do not resterilize
	Caution		Fragile, handle with care
	Upward		Keep dry
LOT	Batch code		Use-by date
	Manufacturer	EU REP	Authorized Representative in the European Union
	Date of manufacture		Consult instructions for use
n	Stacking limit by number		Keep away from sunlight
CE	" CE " Marking	MD	Medical device
UDI	Unique Device Identification	LATEX	Contains latex

【**Storage and transportation**】

1. The product should be stored in a room with a relative humidity of not more than 80%, good ventilation and free of corrosive gases.
2. The product should be protected from heavy pressure, strong vibration, direct sunlight, high temperature, moisture, rain and corrosive gases or liquids during transportation.

【**Sterilization**】 EO sterilization

【**Shelf Life**】2 years after sterilization



Name: Jiangsu HanHeng Medical Technology Co., Ltd.
Add: #118 Xinyuan 3rd Road, Xinbei District, Changzhou, 213031 Jiangsu P.R. China
Tel: 0086-519-85559288
Web: https://hanheng-medical.com

EU REP

Name: Luxus Lebenswelt GmbH
Add: Kochstr. 1, 47877, Willich, Germany
Contact Person: Liu Sun
Tel/Fax: 0049-1715605732
E-mail: Info.m@luxuslw.de

【**Version**】 2025.03.26 B/1

CE 0197