INSTRUCTIONS FOR USE

Umbilical Cord Clamp

【Product Name】 Umbilical Cord Clamp

【Product Model】 HH841001; HH841002; HH841003; HH841004

【Product structure】

Umbilical cord clamp is made of PP (Polypropylene) or polyethylene (PE) material.

【Intended Use】

It is used to clamp over the newborn infant's umbilical cord at delivery.

【Contraindication】

Non

【Steps to use】

1. Wear sterile gloves, open the pouch package and remove the product.
2. Clamp the umbilical cord Clamp tightly into the umbilical cord cutter.
3. After the newborn child is born, the umbilical cord clamp is clamped and cut into the umbilical cord along the umbilical cord clamp of the newborn child.

【Precaution】

1. Please check the pouch package and the product carefully before use. If either of them is damaged, it is strictly forbidden to use.
2. Please read the instruction manual carefully, pay attention to the expiry date of the product before use. Do not use it after the deadline.
3. This product is disposable . It should be strictly in accordance with the aseptic operation specifications. Please destroy it immediately after use. It is forbidden to reuse.
4. The disposal of wastes shall be carried out in accordance with national environmental protection laws and regulations.

【Definitions of Signs】

	Sterilized using ethylene oxide		Do not re-use
	Do not use if package is damaged		Do not resterilize
	Caution		Fragile, handle with care
	Upward		Keep dry

	Batch code		Use-by date
	Manufacturer		Authorized Representative in the European Union
	Date of manufacture		Consult instructions for use
	Temperature limit		Keep away from sunlight
	Stacking limit by number		Medical device
	" CE " Marking		Unique Device Identification

【Storage and transportation】

1. The product should be stored in a relative humidity of no more than 80%, no corrosive gases, dry, well-ventilated room.

2. The product should be protected from heavy pressure, strong vibration, direct sunlight, high temperature, moisture, rain and corrosive gases or liquids during transportation.

【Sterilization】 EO sterilization

【Shelf Life】5 years after sterilization



Name: Jiangsu HanHeng Medical Technology Co., Ltd.

Add: #118 Xinyuan 3rd Road, Xinbei District, Changzhou, 213031 Jiangsu, P.R. China

Tel: 0086-519-85559288

Web: https://hanheng-medical.com

EU REP Name: Luxus Lebenswelt GmbH

Add: Kochstr. 1, 47877, Willich, Germany

Contact Person: Liu Sun

Tel/Fax: 0049-1715605732

E-mail: Info.m@luxuslw.de

【**Version**】2026.05.26 B/1

CE 0197